



GOOD BIOTECH CORP.

Website : www.good-biotech.com

E-mail : technical@good-biotech.com

38, 34th Road Taichung Industrial Park

Taichung 40768 Taiwan

Tel: +886-4-23596873

Fax:+886-4-23590410

ISO9001:2000/ISO13485:2003 & IVD GMP CERTIFIED

SPECIFICATION SHEET

PRODUCT:	Duck IgY(Δ Fc) to alpha-lipase
CATALOG NUMBER	D0702-2B
HOST ANIMAL:	Duck
DESCRIPTION:	IgY(Δ Fc) to alpha-lipase is developed in duck by using purified porcine alpha-lipase as immunogen. The antibody is purified from egg yolk. The product is supplied as a solution in 10mM phosphate buffer, 0.14M NaCl, pH 7.0, containing 0.02% NaN ₃ as a preservative.
IMMUNOGEN:	Alpha-lipase purified from porcine pancreas
CONCENTRATION:	≥ 10 mg/ml, The concentration was determined by 280 nm by spectrophotometer, using 1.4 as the extinction coefficient
PURITY:	Determined by SDS-PAGE (blue stain) , showing that the IgY(Δ Fc) constitute was $\geq 80\%$ of total protein, and the total immunoglobulin IgY+IgY(Δ Fc) constitute was $\geq 90\%$ of total protein.
STORAGE:	For continuous use, store at 0~5°C. For extended storage: -20°C. Repeated freezing and thawing is not recommended.
SPECIFICITY:	This antibody specifically recognized alpha-lipase in western blot analysis and ELISA
TITER:	ELISA titer: $\geq 500x$ ELISA titer is defined as the dilution of antibody (1mg/ml) sufficient to give a change in absorbance of 1.0 at 450nm. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.
APPLICATIONS:	Good to be used in immunoassay method e.g. ELISA, Western Blot assay, Immunoelectrophoresis, Immunodiffusion assay <i>et al.</i>

RELATIVE
PRODUCTS:

For *in vitro* research and the manufacturing of diagnostic only.