



GOOD BIOTECH CORP.

Website : www.good-biotech.com

E-mail : technical@good-biotech.com

38, 34th Road Taichung Industrial Park

Taichung 40768 Taiwan

Tel: +886-4-23596873

Fax:+886-4-23590410

ISO9001:2000/ISO13485:2003 & IVD GMP CERTIFIED

SPECIFICATION SHEET

PRODUCT: Rabbit serum to **Human normal serum**

CATALOG NUMBER R0001-0A

HOST ANIMAL: Rabbit

DESCRIPTION: Antiserum is developed in rabbit by using purified normal human serum as immunogen, containing 0.1% NaN₃ as a preservative.

IMMUNOGEN: Normal human serum

CONCENTRATION: Undiluted, Neat whole serum

STORAGE: For continuous use, store at 0~5°C.
For extended storage: -20°C.
Repeated freezing and thawing is not recommended.

SPECIFICITY: N.D.

TITER: Gel double diffusion titer : **1 : \geq 32**
Gel double diffusion titer is defined as the dilution of antiserum sufficient to form a precipitate arc with 5µg antigen. The antibody determined here was applied 5µl per well. The result was visualized by coomassie blue.

APPLICATIONS: Good to be used in immunoassay method e.g. Immunoturbidimetry, ELISA, Western blot assay, Immunoelectrophoresis, Immunodiffusion assay *et al.*

RELATIVE PRODUCTS: Rabbit serum to Human serum C4 free (**R0002-0A**)

For *in vitro* research and manufacturing of diagnostic only.