



GOOD BIOTECH CORP.

Website : www.good-biotech.com

E-mail : technical@good-biotech.com

38, 34th Road Taichung Industrial Park

Taichung 40768 Taiwan

Tel: +886-4-23596873

Fax:+886-4-23590410

ISO9001:2000/ISO13485:2003 & IVD GMP CERTIFIED

SPECIFICATION SHEET

PRODUCT:	Rabbit antiserum to Human serum albumin
CATALOG NUMBER	R0101-0A
HOST ANIMAL:	Rabbit
DESCRIPTION:	Antiserum is developed in rabbit by using purified human serum albumin as immunogen, containing 0.1% NaN ₃ as a preservative.
IMMUNOGEN:	Purified human serum albumin
CONCENTRATION:	Undiluted, Neat whole serum
STORAGE:	For continuous use, store at 0~5°C. For extended storage: -20°C. Repeated freezing and thawing is not recommended.
SPECIFICITY:	A single arc was shown by Immuno-electrophoresis while the antiserum was reacting with part-pure albumin and normal human serum. Similar specific reactions were observed when the antiserum was tested by western blot assay.
TITER:	Becker titer : ≥ 1.0 mg/ml (antigen/antiserum)
APPLICATIONS:	Good to be used in immunoassay method e.g. Immuno-electrophoresis, Western Blot assay, Immunodiffusion assay, ELISA et al.
RELATIVE PRODUCTS:	Duck IgY(ΔFc) to Human serum albumin (D0101-2B) Duck IgY(ΔFc) to Human serum albumin—affinity purify (D0101-2D) Rabbit IgG to Human serum albumin—affinity purify (R0101-1D) Rabbit IgG to Human serum albumin (R0101-1E) Microalbumin TIA (5010) mAlb Calibrator Set (5511) mAlb Control-L (6011) mAlb Control-H (6012)

For *in vitro* research and manufacturing of diagnostic only.