



**GOOD BIOTECH CORP.**

Website : [www.good-biotech.com](http://www.good-biotech.com)

E-mail : [technical@good-biotech.com](mailto:technical@good-biotech.com)

38, 34<sup>th</sup> Road Taichung Industrial Park

Taichung 40768 Taiwan

Tel: +886-4-23596873

Fax:+886-4-23590410

**ISO9001:2000/ISO13485:2003 & IVD GMP CERTIFIED**

## SPECIFICATION SHEET

PRODUCT:	Rabbit serum to <b>Human transferrin</b>
CATALOG NUMBER	R0105-0A
HOST ANIMAL:	Rabbit
DESCRIPTION:	Antiserum is developed in rabbit by using purified human transferrin as immunogen, containing 0.1% NaN <sub>3</sub> as a preservative.
IMMUNOGEN:	Purified human transferrin
CONCENTRATION:	Undiluted, Neat whole serum
STORAGE:	For continuous use, store at 0~5°C. For extended storage: -20°C. Repeated freezing and thawing is not recommended.
SPECIFICITY:	A single arc was shown while the antiserum was reacting with part-pure human transferrin and normal human serum by Immunoelectrophoresis. Similar specific reactions were observed when the antiserum was tested by western blot assay.
TITER:	Becker titer : $\geq 1.0$ mg/ml ( antigen/antibody )
APPLICATIONS:	Good to be used in immunoassay method e.g. Immunoturbidimetry, ELISA, Western blot assay, Immunoelectrophoresis, Immunodiffusion assay <i>et al.</i>
RELATIVE PRODUCTS:	Duck IgY( $\Delta$ Fc) to Human transferrin ( <b>D0105-2B</b> ) Duck IgY( $\Delta$ Fc) to Human transferrin—affinity purify ( <b>D0105-2D</b> ) Rabbit IgG to Human transferrin—affinity purify ( <b>R0105-1D</b> ) Transferrin TIA ( <b>KT-0105-A</b> ) Transferrin Calibrator Set ( <b>CA-0105-A</b> ) Transferrin Control-L ( <b>CL-0105-A</b> ) Transferrin Control-M ( <b>CL-0105-B</b> ) Transferrin Control-H ( <b>CL-0105-C</b> )

For *in vitro* research and manufacturing of diagnostic only.