



GOOD BIOTECH CORP.

Website : www.good-biotech.com

E-mail : technical@good-biotech.com

38, 34th Road Taichung Industrial Park

Taichung 40768 Taiwan

Tel: +886-4-23596873

Fax:+886-4-23590410

ISO9001:2000/ISO13485:2003 & IVD GMP CERTIFIED

SPECIFICATION SHEET

PRODUCT:	Rabbit serum to Human hepsin
CATALOG NUMBER	R0203-0A
HOST ANIMAL:	Rabbit
DESCRIPTION:	Antiserum is developed in rabbit by using synthetic peptide as immunogen, containing 0.1% NaN ₃ as a preservative.
IMMUNOGEN:	synthetic peptide
CONCENTRATION:	Undiluted, Neat whole serum
STORAGE:	For continuous use, store at 0~5°C. For extended storage: -20°C. Repeated freezing and thawing is not recommended.
SPECIFICITY:	This antiserum specifically recognized recombinant human hepsin in western blot analysis; cross-reactivity to hepsin of other species has not yet been tested.
TITER:	ELISA titer : ≥50000x ELISA titer is defined as the dilution of antiserum sufficient to give a change in absorbance of 1.0 at 450nm. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.
APPLICATIONS:	Good to be used in immunoassay method e.g. Immunohistochemistry, Cell flowcytometry, ELISA, Western Blot assay <i>et al.</i> Other applications for use of this antibody are available, but the performances should be estimated by users.
RELATIVE PRODUCTS:	Duck IgY(ΔFc) to Human hepsin—affinity purified (D0203-2D) Rabbit IgG to Human hepsin—affinity purified (R0203-1D)

For *in vitro* research and manufacturing of diagnostic only.