



GOOD BIOTECH CORP.

Website : www.good-biotech.com

E-mail : technical@good-biotech.com

38, 34th Road Taichung Industrial Park

Taichung 40768 Taiwan

Tel: +886-4-23596873

Fax: +886-4-23590410

ISO9001:2000/ISO13485:2003 & IVD GMP CERTIFIED

SPECIFICATION SHEET

| | |
|--------------------|---|
| PRODUCT: | Rabbit serum to Human C-reactive protein |
| CATALOG NUMBER | R0301-0A |
| HOST ANIMAL: | Rabbit |
| DESCRIPTION: | Antiserum is developed in rabbit by using purified human CRP as immunogen, containing 0.1% NaN ₃ as a preservative. |
| IMMUNOGEN: | Purified human C-reactive protein |
| CONCENTRATION: | Undiluted, Neat whole serum |
| STORAGE: | For continuous use, store at 0~5°C. For extended storage: -20°C. Repeated freezing and thawing is not recommended. |
| SPECIFICITY: | A single arc was shown while the antiserum was reacting with part-pure human C-reactive protein by Immunoelectrophoresis. Nothing was shown when the antiserum reacted with normal human serum. Similar specific reactions were observed when the antiserum was tested by western blot assay. |
| TITER: | Becker titer : ≥ 2.0 mg/ml (antigen/antibody) |
| APPLICATIONS: | Good to be used in immunoassay method e.g. Immunoturbidimetry, ELISA, Western blot assay, Immunoelectrophoresis, Immunodiffusion assay <i>et al.</i> |
| RELATIVE PRODUCTS: | Duck IgY(Δ Fc) to Human C-reactive protein (D0301-2B) Goat serum to Human C-reactive protein—gamma fraction (G0301-1B) Rabbit IgG to Human C-reactive protein (R0301-1E) CRPex-HS CRP LIT Assay (1010) CRPex-HS CRP Calibrator Set (1510) CRPex-BR CRP LIT Assay (2010) CRPex-BR CRP Calibrator Set (2510) CRPex CRP Control-L (3011) CRPex CRP Control-M (3012) CRPex CRP Control-H (3013) Duet hs-CRP LIT Assay (4010) Duet hs-CRP Calibrator Set (4511) Triplex CRP LIT Assay (KT-0301-D) |

For *in vitro* research and manufacturing of diagnostic only.