



GOOD BIOTECH CORP.

Website : www.good-biotech.com

E-mail : technical@good-biotech.com

38, 34th Road Taichung Industrial Park

Taichung 40768 Taiwan

Tel: +886-4-23596873

Fax:+886-4-23590410

ISO9001:2000/ISO13485:2003 & IVD GMP CERTIFIED

SPECIFICATION SHEET

PRODUCT:	Rabbit IgG to duck IgY(ΔFc) —HRP conjugated
CATALOG NUMBER	R8102-11
HOST ANIMAL:	Rabbit
DESCRIPTION:	Affinity purified rabbit IgG to duck IgY(Δ Fc) is purified by duck IgY(Δ Fc) coupling column. The product was conjugated with horse radish peroxidase. The product is supplied as a solution in 50mM Tris, 0.14M sodium chloride, pH 7.9.
IMMUNOGEN:	Purified duck IgY(Δ Fc)
STORAGE:	Please store at -20°C
SPECIFICITY:	This antibody specifically recognized IgY(Δ Fc) in immunoelectrophoresis and western blot assay.
TITER:	ELISA titer : $\geq 10000x$ ELISA titer is defined as the dilution of antibody sufficient to give a change in absorbance of 1.0 at 450nm. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.
APPLICATIONS:	This antibody is suitable as a 2 nd antibody while being 5000x-diluted for the usage in ELISA and western blot assay.
RELATIVE PRODUCTS:	Rabbit antiserum to Duck Immunoglobulin Y(Δ Fc) (R8102-0A) Rabbit IgG to Duck Immunoglobulin Y(Δ Fc)--Affinity purified (R8102-1D) Rabbit IgG to Duck Immunoglobulin Y(Δ Fc) (R8102-1E)

For *in vitro* research and manufacturing of diagnostic only.