



**GOOD BIOTECH CORP.**

Website : [www.good-biotech.com](http://www.good-biotech.com)

E-mail : [technical@good-biotech.com](mailto:technical@good-biotech.com)

38, 34<sup>th</sup> Road Taichung Industrial Park

Taichung 40768 Taiwan

Tel: +886-4-23596873

Fax:+886-4-23590410

**ISO9001:2008/ISO13485:2012 & IVD GMP CERTIFIED**

---

## SPECIFICATION SHEET

PRODUCT:	Goat antiserum to <b>Human C-reactive protein</b>
CATALOG NUMBER	G0301-0A
LOT NUMBER:	16GS102
HOST ANIMAL:	Goat
DESCRIPTION:	Antiserum is developed in goat by using purified human CRP as immunogen, containing 0.1% NaN <sub>3</sub> as a preservative.
IMMUNOGEN:	Purified human C-reactive protein
CONCENTRATION:	Undiluted, Neat whole serum
STORAGE:	For continuous use in a short period, store at 0~5°C. For extended storage: -20°C. Repeated freezing and thawing is not recommended.
SPECIFICITY:	A single arc was shown while the antiserum was reacting with part-pure human C-reactive protein by Immuno-electrophoresis. Nothing was shown when the antiserum reacted with normal human serum. Similar specific reactions were observed when the antiserum was tested by western blot assay.
TITER:	Becker titer : <b>3.5</b> mg/ml ( antigen/antibody )
SHELF LIFE	3 years after dispatch
APPLICATIONS:	Good to be used in immunoassay method e.g. Immunoturbidimetry, ELISA, Western blot assay, Immuno-electrophoresis, Immunodiffusion assay <i>et al.</i>

**For *in vitro* research and manufacturing of diagnostic only.**

---